

## ABSTRACT OF THE DISCLOSURE

Applicable to humans, animals and fish, a method for injecting thoroughly diffused ambient  
5 air or disinfectant into water prior to its delivery into a therapy tank plus an underwater PZT  
probe transmission of separate stable and transient cavitation signals from which a  
microcomputer determines, 1) the average number of transducer generated sinoidal equal  
amplitude alternating compression and rarefaction ultrasonic acoustic pressure waves cycles  
necessary to create inertial and/or transient cavitation and the required number of rectified  
10 sinoidal equal amplitude ultrasonic compression acoustic pressure waves necessary to  
suppress the inertial/transient cavitation and thereby maintain stable cavitation for cleaning  
and open-wound therapy treatment for 15-minutes, (or greater) time periods and, 2) the  
necessary dilution of water and disinfectant and its activation by dual-mode transient  
cavitation to kill the pathogens shed by the "patient" following "patient" cleaning or wound-  
15 therapy treatment.